



Development and validation of a new homogeneous immunoassay for the detection of UR-144 metabolites in urine

Guohong Wang*, Kim Huynh, Rehka Barhate, Tony Cam, Jialin Liu, and Warren Rodrigues
Immunalysis Corporation, Pomona, CA USA



Introduction and Objective

Synthetic cannabinoids are consistently being synthesized worldwide with the intention of evading detection in biological matrices. The objective of this project was to develop and validate a new high throughput homogeneous enzyme immunoassay (HEIA) for the rapid detection of the urinary metabolites of the latest synthetic cannabinoids such as UR-144 and XLR-11. Despite the fact that UR-144, XLR-11, JWH-018, and AM2201 belong to a drug class that share a core indole-ring moiety, it is challenging to develop an immunoassay to cross react with all of them due to the significant difference in structures of the substituted groups. The current commercially available homogeneous immunoassay targeted at JWH- metabolites has very low cross reactivity for UR-144 and XLR-11, thus it was necessary to develop a new screening method to detect UR-144, XLR-11 and their metabolites in urine.

Advantages

1. Ready to use reagents suitable for high throughput instruments
2. Assay working range : 0 to 40ng/mL, at a cutoff of 10ng/mL
3. Specific assay with accuracy >98% based on 65 urine specimens

Results and Discussion

Precision: Daily Calibration Required

The qualitative precision was determined by assaying calibrators and controls in synthetic urine for 20 days, 2 runs per day in duplicates (N=80). The results are summarized below.

Conc. ng/mL	Mean Conc. (ng/ml)	C.V.%
5 ng/mL (Control LOW)	537.2	4.6
10 ng/mL Calibrator	564.2	5.0
15ng/mL (Control HIGH)	585.8	4.7

Cross Reactivity:

Structurally related compounds that are potentially found in urine were tested using the 10 ng/mL cutoff calibrator.

Structurally Related Compounds			
Compound	Analyte Conc. (ng/mL)	UR-144 N-pentanoic acid (ng/mL)	Cross-Reactivity (%)
UR-144 N-pentanoic acid	10	10	100
UR-144	20	10	50
UR-144 N-(5hydroxypentyl) β -d-glucuronide	30	10	33
UR-144 N-(5-bromopentyl)	25	10	40
UR-144 N-(5-chloropentyl)	20	10	50
UR-144 N-heptyl	40	10	25
UR-144 N-(5hydroxypentyl) metabolite	20	10	50
XLR-11	20	10	50
XLR-11 N-(4-pentenyl)	20	10	50
XLR-11 N-(4-hydroxypentyl) metabolite	70	10	14

Structurally Related Compounds

Compound	Analyte Conc. (ng/mL)	UR-144 N-pentanoic acid (ng/mL)	Cross-Reactivity (%)
AB-005	30	10	33
A-834735	20	10	50
JWH-250 5 hydroxyindole metabolite	100,000	10	ND
JWH-250-N(hydroxypentyl) metabolite	20,000	10	0.05
JWH-250-N(4-hydroxypentyl) metabolite	50,000	10	0.02
JWH-250-N(5-carboxypentyl) metabolite	50,000	10	0.02
RCS-4-2 methoxy isomer	10,000	10	0.10
JWH-250	100,000	10	ND
JWH-210	100,000	10	ND
JWH-018 N(5-hydroxypentyl)metabolite	3,000	10	0.3
AM2201-6-hydroxyindole metabolite	100,000	10	ND
JWH-007	100,000	10	ND
JWH-019	100,000	10	ND
AM-2232	100,000	10	ND
JWH-081	100,000	10	ND
AM-2201-N-4-hydroxypentyl metabolite	100,000	10	ND
AM-2233	10,000	10	0.1%
1-Naphthoyl indole	100,000	10	ND
JWH-073 6-hydroxyindole metabolite	100,000	10	ND
JWH-073	100,000	10	ND
JWH-122	100,000	10	ND
JWH-018 N(5-hydroxypentyl β -d-glucuronide	100,000	10	ND
JWH-201	100,000	10	ND
AM-2201	100,000	10	ND
Cannabipiperidiethanone	50,000	10	0.02%
JWH-018 4-hydroxyindole	100,000	10	ND
JWH-022	100,000	10	ND
JWH-073 N-butanoic acid	100,000	10	ND
JWH-018 5-hydroxyindole	100,000	10	ND
3-1 naphthoyl-1H-indole	100,000	10	ND
JWH-018 pentanoic acid	100,000	10	ND
JWH-018	100,000	10	ND

Authentic specimens:

- 65 urine specimens previously confirmed by an outside laboratory by LC-MS/MS were analyzed with this Immunalysis UR-144 PA EIA assay
- Cutoff concentration: 10 ng/mL for both EIA and confirmation method
- 40 specimens were negative by both methods
- 24 specimens were positive by both methods
- The sensitivity, specificity and accuracy were 96%, 100%, and 98%, respectively

		LC-MS/MS	
		+	-
HEIA	+	24	0
	-	1	40

- Negative result: absorbance rate reading just below cutoff

Summary

A high throughput HEIA has been developed for the detection of UR-144 metabolites in human urine which correlates well with LC-MS/MS.

References

- Grigoryev, A., Kavanagh, P., Melnik, P., Savchuk, S., Simonov, A.; Gas and Liquid Chromatography-Mass Spectrometry Detection of the Urinary Metabolites of UR-144 and Its Major Pyrolysis Product; *Journal of Analytical Toxicology*, April 16, 2013: doi:10.1093/jat/bkt028.

SOFT, Orlando, FL 2013