



User Friendly – Flexible – Reliable

The **ImmTox™ Benchtop Chemistry Analyzer** offers a comprehensive menu of toxicology assays. **ImmTox** is an ideal solution for physician offices and treatment centers needing timely, relevant, lab quality results.

- **Throughput:** 240 *dual* reagent tests per hour
- **Patience Sample Capacity:** 40 on board
- **Reagent Capacity:** 24 assays on board
- **Bi-directional LIS communication**
- **Fully traceable** results and quality control data
- **Continuous access** to samples and reagents without interruption of the testing process
- **Reusable cuvettes** – more cost effective

LIS Software

Developed Specifically for the Physician's Office Laboratory and Treatment Centers



- Intuitive user interface
- Easy management of patient records
- Flexibility to generate summary reports for patients and quality control
- EMR integration options

COMPREHENSIVE TEST MENU



Drug/Analyte	Brand Name or Street Name
Opiates and Opioids	
Buprenorphine	Suboxone®, Subutex®, Temgesic®
Codeine and Morphine	Hycodan®, Lorcet®, Avinza®, Roxanol®, MS Contin®
Fentanyl	Actiq®, Duragesic®, Onsolis®
Meperidine	Demerol®, Pethidine
Methadone	Dolophine®, Methadose®
Oxycodone and Hydrocodone	OxyContin®, Percocet®
Tapentadol	Nucynta®
Tramadol	Ultram®, Ryzolt®
Benzodiazepines	
	Xanax®, Valium®, Ativan®, Serax®
Barbiturates	
	Butisol®, Fiorinal®, Phenobarbital
Stimulants	
Amphetamine	Adderall®, Dexedrine®
Methamphetamine	Desoxyn®
Nicotine (Cotinine)	Nicoderm®, Nicorette®
Muscle Relaxants	
Carisoprodol and Meprobamate	Soma®, Miltown®, Equanil®
Sedative/Hypnotics	
Zolpidem	Ambien®, Stilnox®, Sublinox®
Tricyclic Antidepressants	
	Tofranil®, Elavil®, Anafranil®
Illicit Drugs	
Cocaine (Benzoylecgonine)	Coke, Blow, Dream, Snow
MDMA (Methylenedioxymethamphetamine)	Adam, Ecstasy, XTC
MDA (Methylenedioxyamphetamine)	Eve, Love Drug
Marijuana (THC metabolite)	Pot, Weed, Grass, Dope, Reefer, Ganja
Phencyclidine (PCP)	Angel Dust
Alcohol / Alcohol Metabolites	
Ethyl Alcohol	
Ethyl Glucuronide (EtG)	

Site Requirements: **Water Consumption:** 3.5 liters (<1 gallon) per hour – no need for external water or drainage connections
Electrical: Standard 110V outlet • **Dimensions:** 32”(W) x 26”(D) x 21”(H) – 38”(H) with open cover



www.immunalysis.com

Tel 909.482.0840 | Toll Free 888.664.8378 | Fax 909.482.0850

©2014 Immunalysis Corporation. All rights reserved. Immunalysis, Immunalysis logo and ImmTox are trademarks of Immunalysis Corporation. MKT-1070 REV. A