




BENCHTOP TOXICOLOGY ANALYZER

# IMMTOX™ 270

Small, simple  
and smart

The ImmTox™ 270 benchtop analyzer offers comprehensive toxicology screening solutions for physician offices, treatment centers and independent laboratories.

- Compact design with easy-to-use interface
- Comprehensive test menu CLIA<sup>1</sup> categorized as moderate complexity
- Low cutoff concentrations
- Fast results
- High throughput
- Walkaway operating functionality
- Minimal maintenance



## FEATURES AND BENEFITS

- 270 dual reagent tests per hour
- Reagent capacity: 24 assays
- Onboard patient sample capacity: up to 30
- Ideally sized for 200-1200 patient samples per month
- Intuitive, touch-screen user interface
- Bi-directional LIS communication
- Fully traceable results and quality control data
- Continuous access to samples without interruption of the testing process
- Reusable cuvettes and wedges
- Fast start-up
- Strong technical support

## LIS COMMUNICATION SOFTWARE AVAILABLE

- Easy management of patient records
- Flexibility to generate summary reports for patients and quality controls
- EMR<sup>2</sup> integration options
- HL7<sup>3</sup> interface

## SPECIFICATIONS

- Dimensions: 32”(W) x 27”(D) x 21”(H) – 38”(H) with open cover
- Electrical: Standard 110V outlet
- Self-contained analyzer
  - Water consumption: 3.8 liters
  - External water or drainage connections not required
  - Easy to install and operate

## TEST MENU

URINE DRUG SCREENING REAGENTS	CUTOFF LEVEL (ng/mL)	MODERATE COMPLEXITY
6-acetylmorphine (6-AM Heroin metabolite)	10	☑
Amphetamine	500 or 1,000	☑
Barbiturates	200	☑
Benzodiazepines	200	☑
Benzoylcegonine (Cocaine metabolite)	150 or 300	☑
Buprenorphine	5	☑
Cannabinoids (THC)	50	☑
Carisoprodol	100	
EDDP (Methadone metabolite)	100 or 300	☑
Ethyl Alcohol	100 mg/dL	
Ethyl Glucuronide (EtG)	500 or 1,000	
Fentanyl (SEFRIA™)	1	☑
Fentanyl (HEIA™)	2	
Ketamine	100	
Meperidine	200	
Methamphetamine	500 or 1,000	☑
Opiates	300 or 2,000	☑
Oxycodone	100 or 300	☑
Phencyclidine (PCP)	25	☑
Synthetic Cannabinoids-1 (JWH-018, JWH-073, AM-2201 and its metabolites)	10	
Synthetic Cannabinoids-2 (UR-144, XLR11 and its metabolites)	10	
Synthetic Cannabinoids-3 (AB-PINACA, ADB-PINACA and its metabolites)	10	
Tapentadol	200	
Tramadol	200	☑
Zolpidem	20	

1. Clinical Laboratory Improvement Amendments (CLIA) 2. Electronic Medical Records (EMR) 3. Health Level Seven International, Inc. (HL7®)

855-425-9428 | [CLS\\_SALES@ABBOTT.COM](mailto:CLS_SALES@ABBOTT.COM)

© 2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. SEFRIA is a trademark of Advanced ImmunoAssays, LLC. Any photos displayed are for illustrative purposes only. COL-02948 v2 01/22

